



#5

HEGAR**DILATOR**

Length: 15/16mm, 8in
Author: Alfred Hegar

EST.
1879



Dilator Hegar Uterine 15/16mm, 8in

Common Uses:

Penile implant procedure and Cervical procedure

Properties:

- Common tip diameter - 16mm to 26mm
- Common Length - 7in to 13in

Variations:

- Single ended
- Double ended

Function:

The dilator is used for assessing the length and planes of the uterus, its direction and often it is used for inducing dilatation of the uterine cavity.

Did you know?

Hegar dilators are used to treat vaginismus, also known as genito-pelvic pain disorder.

